

AP0832

Leader in Biomolecular Solutions for Life Science



Phospho-eEF2-T56 Rabbit pAb

Catalog No.: AP0832

Basic Information

Observed MW

95kDa

Calculated MW

95kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

This gene encodes a member of the GTP-binding translation elongation factor family. This protein is an essential factor for protein synthesis. It promotes the GTP-dependent translocation of the nascent protein chain from the A-site to the P-site of the ribosome. This protein is completely inactivated by EF-2 kinase phosphorylation.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

1938

Swiss Prot

P13639

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

EF2; EF-2; EEF-2; SCA26; Phospho-eEF2-T56

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from 293T cells, using Phospho-eEF2-T56 Rabbit pAb (AP0832) at 1:3000 dilution. 293T cells were treated with 20% FBS at 37°C for 20 minutes after serum-starvation overnight

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 5s.