

AP0926

Leader in Biomolecular Solutions for Life Science



Phospho-ARHGEF2-S886 Rabbit pAb

Catalog No.: AP0926

Basic Information

Observed MW

120kDa

Calculated MW

112kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Rho GTPases play a fundamental role in numerous cellular processes that are initiated by extracellular stimuli that work through G protein coupled receptors. The encoded protein may form complex with G proteins and stimulate rho-dependent signals. Alternatively spliced transcript variants encoding different isoforms have been identified.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9181

Swiss Prot

Q92974

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

GEF; Lfc; P40; GEFH1; LFP40; GEF-H1; NEDMHM; Phospho-ARHGEF2-S886

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from C6 cells, using Phospho-ARHGEF2-S886 Rabbit pAb (AP0926) at 1:1000 dilution. C6 were treated by CIP (20ul CIP for each 400ul cell lysate) at 37°C for 1 hour.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25μg per lane.

Blocking buffer: 3% BSA.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.