

AP1237

Leader in Biomolecular Solutions for Life Science



Phospho-IKK β -Y199 Rabbit pAb

Catalog No.: AP1237

Basic Information

Observed MW

87kDa

Calculated MW

87kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

The protein encoded by this gene phosphorylates the inhibitor in the inhibitor/NF-kappa-B complex, causing dissociation of the inhibitor and activation of NF-kappa-B. The encoded protein itself is found in a complex of proteins. Several transcript variants, some protein-coding and some not, have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

3551

Swiss Prot

O14920

Immunogen

A synthetic phosphorylated peptide around Y199 of human IKK β (NP_001547.1).

Synonyms

IKK2; IKKB; IMD15; IMD15A; IMD15B; NFKB1KB; IKK-beta; Phospho-IKK β -Y199

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.