

RK05632

Leader in Biomolecular Solutions for Life Science



[One Step] CCT8 Antibody Kit

Catalog No.: RK05632

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

60kDa

Calculated MW

60kDa

Category

Antibody kit

Background

This gene encodes the theta subunit of the CCT chaperonin, which is abundant in the eukaryotic cytosol and may be involved in the transport and assembly of newly synthesized proteins. Alternative splicing results in multiple transcript variants of this gene. A pseudogene related to this gene is located on chromosome 1.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05632-1	CCT8 Rabbit pAb	60kDa
RK05632-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	60kDa

Immunogen Information

Gene ID

10694

Swiss Prot

P50990

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 248-497 of human CCT8 (NP_001269837.1).

Synonyms

Cctq; PRED71; D21S246; C21orf112; CCT8