

RK05645

Leader in Biomolecular Solutions for Life Science



[One Step] eIF5A Antibody Kit

Catalog No.: RK05645

Basic Information

Applications

Cross-Reactivity

Human, Mouse

Observed MW

17kDa

Calculated MW

17kDa

Category

Antibody kit

Background

Enables U6 snRNA binding activity and protein N-terminus binding activity. Involved in several processes, including cellular response to virus; positive regulation of intrinsic apoptotic signaling pathway by p53 class mediator; and tumor necrosis factor-mediated signaling pathway. Located in annulate lamellae; cytoplasm; and nucleus. Part of nuclear pore.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05645-1	eIF5A Rabbit pAb	17kDa
RK05645-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	17kDa

Immunogen Information

Gene ID

1984

Swiss Prot

P63241

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-154 of human eIF5A (NP_001137234.1).

Synonyms

FABAS; EIF-5A; EIF5A1; eIF-4D; eIF5A1; eIF5A