

RK05663

Leader in Biomolecular Solutions for Life Science



[One Step] LSM4 Antibody Kit

Catalog No.: RK05663

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

16kDa

Calculated MW

15kDa

Category

Antibody kit

Background

This gene encodes a member of the LSm family of RNA-binding proteins. LSm proteins form stable heteromers that bind specifically to the 3'-terminal oligo(U) tract of U6 snRNA and may play a role in pre-mRNA splicing by mediating U4/U6 snRNP formation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05663-1	LSM4 Rabbit pAb	15kDa
RK05663-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	16kDa

Immunogen Information

Gene ID

25804

Swiss Prot

Q9Y4Z0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-139 of human LSM4 (NP_036453.1).

Synonyms

GRP; YER112W; LSM4