

RK05669

Leader in Biomolecular Solutions for Life Science



## [One Step] NDUFA7 Antibody Kit

Catalog No.: RK05669

### Basic Information

#### Applications

#### Cross-Reactivity

Human, Mouse, Rat

#### Observed MW

13kDa

#### Calculated MW

13kDa

#### Category

Antibody kit

### Background

This gene encodes a subunit of NADH:ubiquinone oxidoreductase (complex I), which is a multiprotein complex located in the inner mitochondrial membrane. Complex I functions in the transfer of electrons from NADH to the respiratory chain.

### Product Information

#### Source

Rabbit

#### Purification

#### Storage

Avoid repeated freeze-thaw cycles.

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05669-1	[KO Validated] NDUFA7 Rabbit pAb	13kDa
RK05669-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	13kDa

### Immunogen Information

#### Gene ID

4701

#### Swiss Prot

O95182

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human NDUFA7 (NP\_004992.2).

#### Synonyms

B14.5a; CI-B14.5a; A7