

RK05684

Leader in Biomolecular Solutions for Life Science



[One Step] PQBP1 Antibody Kit

Catalog No.: RK05684

Basic Information

Applications

Cross-Reactivity

Human, Mouse, Rat

Observed MW

30kDa

Calculated MW

30kDa

Category

Antibody kit

Background

This gene encodes a nuclear polyglutamine-binding protein that is involved with transcription activation. The encoded protein contains a WW domain. Mutations in this gene have been found in patients with Renpenning syndrome 1 and other syndromes with X-linked cognitive disability. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

| Catalog No. | Product Name | Dilutions |
|-------------|---|-----------|
| RK05684-1 | PQBP1 Rabbit pAb | 30kDa |
| RK05684-2 | HRP-conjugated Goat anti-Rabbit IgG (H+L) | 30kDa |

Immunogen Information

Gene ID

10084

Swiss Prot

O60828

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human PQBP1 (NP_005701.1).

Synonyms

SHS; MRX2; MRX55; MRX53; MRX58; NPW38; RENS1; PQBP1