

RK05715

Leader in Biomolecular Solutions for Life Science



[One Step] TIMM17A Antibody Kit

Catalog No.: RK05715

Basic Information

Applications

Cross-Reactivity

Human, Mouse

Observed MW

18kDa

Calculated MW

18kDa

Category

Antibody kit

Background

Predicted to contribute to protein transmembrane transporter activity. Predicted to be involved in protein import into mitochondrial matrix. Located in mitochondrial inner membrane and nucleoplasm. Is integral component of mitochondrial inner membrane. Part of TIM23 mitochondrial import inner membrane translocase complex.

Product Information

Source

Rabbit

Purification

Storage

Avoid repeated freeze-thaw cycles.

Contact



www.abclonal.com

Component & Recommended Dilutions

Catalog No.	Product Name	Dilutions
RK05715-1	TIMM17A Rabbit pAb	18kDa
RK05715-2	HRP-conjugated Goat anti-Rabbit IgG (H+L)	18kDa

Immunogen Information

Gene ID

10440

Swiss Prot

Q99595

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-171 of human TIMM17A (NP_006326.1).

Synonyms

TIM17; TIM17A; TIMM17A