

RMK0240

Leader in Biomolecular Solutions for Life Science



# Rabbit anti-Human sE-Selectin/CD62E mAb(CAP)

Catalog No.: RMK0240

## Basic Information

### Catagory

Matched Antibody Pair

### Applications

multiplex assay

## Product Information

### Ig Type

Rabbit IgG

### Purification

Affinity purification

### Endotoxin Level

### Storage

Store unopened kit at 2-8 °C for 1 week. If more than one week, please keep the components of the kit according to the instructions. Do not use past kit expiration date.

**Avoid repeated freeze-thaw cycles.**

### Formulation

Supplied as a 0.2um filtered solution in PBS, PH 7.4 containing Human sE-Selectin, CD62E Antibody.

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Background

E-Selectin, also known as Endothelial Leukocyte Adhesion Molecule 1 (ELAM-1), and Leukocyte-Endothelial Cell Adhesion Molecule 2 (LECAM-2), has been designated CD62E. It is a cell-surface glycoprotein expressed by activated endothelium and it mediates the adhesion of blood neutrophils.

## Immunogen Information

### Immunogen

Recombinant Human sE-Selectin/CD62E Protein

## Cross-Reactivity


/

## Assay Applications

Human sE-Selectin/CD62E multiplex assay

	Recommended Concentration	Sample
Multiplex Capture	3-20ug/mL	Rabbit anti-Human sE-Selectin/CD62E mAb(Cat. No.RMK0240)
Multiplex Detection	0.017-2ug/mL	Rabbit anti-Human sE-Selectin/CD62E mAb(Cat. No.RMK0241)
Standard	6.86-5000pg/mL	Recombinant Human sE-Selectin/CD62E Protein

# Validation Data



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.