

RMK0531

Leader in Biomolecular Solutions for Life Science



# Rabbit anti-Human Phospho-Tau-T217 mAb

Catalog No.: RMK0531

## Basic Information

### Catagory

Matched Antibody Pair

### Applications

CLIA

## Product Information

### Ig Type

Rabbit IgG

### Purification

Affinity purification

### Endotoxin Level

### Storage

Store at -20°C.

Avoid repeated freeze-thaw cycles.

### Formulation

Supplied as a 0.2µm filtered solution in PBS with 0.05% Proclin300,PH 7.4.

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Background

## Immunogen Information

### Immunogen

Recombinant Human Phospho-Tau-T217 Protein

## Cross-Reactivity

/


## Assay Applications

Human Phospho-Tau-T217 Sandwich Immunoassay

		Recommended Concentration	Sample
CLIA	Capture	1-4ug/mL	Rabbit anti-Human Phospho-Tau-T217 mAb(Cat. No.RMK0531)
	Detection	0.5-2ug/mL	Rabbit anti-Human Tau mAb(Cat. No.RMK0456)
	Standard	1-5000pg/mL	Recombinant Human Phospho-Tau-T217 Protein

## Validation Data

---



This standard curve is only for demonstration purposes. A standard curve should be generated for each assay.