

RP00070

Leader in Biomolecular Solutions for Life Science



Recombinant Human TNFRSF14/HVEM/CD270 Protein

Catalog No.: RP00070

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	8764	Q92956-1

Tags

C-hFc&His

Synonyms

TNFRSF14; ATAR; CD270; HVEA; HVEM; LI
GHTR; TR2

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Background

Herpesvirus entry mediator (HVEM), also known as tumor necrosis factor receptor superfamily member 14 (TNFRSF14), is a human cell surface receptor of the TNF-receptor superfamily. Two TNF superfamily ligands lymphotoxin α (TNF-β) and LIGHT (TNFSF14) are identified as cellular ligands for HVEM and initiate the positive signaling.

Basic Information

Description

Recombinant Human TNFRSF14/HVEM/CD270 Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Leu 39 - Val 202) of human HVEM (Accession #NP_003811.2) fused with an Fc, 6×His tag at the C-terminus.

Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Recombinant human HVEM at 5 μg/mL (100 μL/well) can bind Biotinylated Recombinant human BTLA with a linear range of 1.5-6 μg/mL.

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.

Reconstitution


Centrifuge the vial before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact




www.abclonal.com

Validation Data



Recombinant Human TNFRSF14/HVEM/CD270 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Recombinant human HVEM at 5 µg/mL (100 µL/well) can bind Biotinylated Recombinant human BTLA with a linear range of 1.5-6 µg/mL.