

RP00094

Leader in Biomolecular Solutions for Life Science



Recombinant Human CSF-2/GM-CSF Protein

Catalog No.: RP00094

Recombinant

2 Publications

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	1437	P04141

Tags

C-His

Synonyms

CSF2;GMCSF

Product Information

Source

HEK293 cells

Purification

≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of granulocytes and macrophages. The active form of the protein is found extracellularly as a homodimer. This gene has been localized to a cluster of related genes at chromosome region 5q31, which is known to be associated with interstitial deletions in the 5q- syndrome and acute myelogenous leukemia. Other genes in the cluster include those encoding interleukins 4, 5, and 13.

Basic Information

Description

Recombinant Human CSF-2/GM-CSF Protein is produced by HEK293 expression system. The target protein is expressed with sequence (Ala18-Glu144) of human GM-CSF/CSF2 (Accession #NP_000749.2) fused with a 6×His tag at the C-terminus.

Bio-Activity

Measured in a cell proliferation assay using TF-1 human erythroleukemic cells. The ED₅₀ for this effect is 30-140 pg/mL, corresponding to a specific activity of 7.14×10^6 - 3.33×10^7 units/mg.

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human CSF-2/GM-CSF Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Recombinant Human CSF-2 stimulates cell proliferation of the TFR1 human erythroleukemic cells. The ED₅₀ for this effect is 30-140 pg/mL, corresponding to a specific activity of 7.14×10^6 - 3.33×10^7 units/mg.