

RP00131

Leader in Biomolecular Solutions for Life Science



# Recombinant Human CEACAM3/CD66d Protein

Catalog No.: RP00131 **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	1084	P40198

### Tags

C-His

### Synonyms

CEACAM3;CD66D;CEA;CGM1;W264;W282

## Product Information

Source	Purification
	≥ 90 % as determined by SDS-PAGE.

### Endotoxin

< 0.1 EU/μg of the protein by LAL method.

### Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

### Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Background

The protein encodes a member of the family of carcinoembryonic antigen-related cell adhesion molecules (CEACAMs), which are used by several bacterial pathogens to bind and invade host cells. The encoded transmembrane protein directs phagocytosis of several bacterial species that is dependent on the small GTPase Rac. It is thought to serve an important role in controlling human-specific pathogens by the innate immune system. Alternatively spliced transcript variants have been described.

## Basic Information

### Description

Recombinant Human CEACAM3/CD66d Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Gly155) of human CEACAM3/CD66d (Accession #NP\_001806.2) fused with a 6×His tag at the C-terminus.

### Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Human CEACAM3 at 1μg/mL (100 μL/well) can bind CEACAM3 Mouse mAb with a linear range of 0.1-2.7 ng/mL.

### Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

### Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.


### Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.


\* For your safety and health, please wear a lab coat and disposable gloves when handling.

## Validation Data

---



Recombinant Human CEACAM3/CD66d Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Human CEACAM3 at 1 $\mu$ g/mL (100  $\mu$ L/well) can bind CEACAM3 Mouse mAb with a linear range of 0.1-2.7 ng/mL.