

RP00152

Leader in Biomolecular Solutions for Life Science



Recombinant Human Low-density lipoprotein receptor/LDLR Protein

Catalog No.: RP00152

Recombinant

2 Publications

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	3949	P01130

Tags

C-His

Synonyms

LDLR; FH; FHC; LDLCQ2; low-density lipoprotein receptor; FH; FHC; LDLCQ2

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 90 % as determined by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Background

The low density lipoprotein receptor (LDLR) family consists of cell surface proteins involved in receptor-mediated endocytosis of specific ligands. Low density lipoprotein (LDL) is normally bound at the cell membrane and taken into the cell ending up in lysosomes where the protein is degraded and the cholesterol is made available for repression of microsomal enzyme 3-hydroxy-3-methylglutaryl coenzyme A (HMG CoA) reductase, the rate-limiting step in cholesterol synthesis. At the same time, a reciprocal stimulation of cholesterol ester synthesis takes place. Mutations in LDLR gene lead to the autosomal dominant disorder, familial hypercholesterolemia.

Basic Information

Description

Recombinant Human Low-density lipoprotein receptor/LDLR Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Ala22-Arg788) of human LDLR/LDL Receptor (Accession #NP_000518.1) fused with a 6×His tag at the C-terminus.

Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Human LDLR at 5 μg/mL (100 μL/well) can bind Human PCSK9 with a linear range of 0.5-294.7 ng/mL.

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.

Reconstitution


Centrifuge the vial before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact




www.abclonal.com

Validation Data



Recombinant Human Low-density lipoprotein receptor/LDLR Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Human LDLR at 5μg/mL (100 μL/well) can bind Human PCSK9 with a linear range of 0.5-294.7ng/mL.