

RP00202

Leader in Biomolecular Solutions for Life Science



Recombinant Human Cripto-1/TDGF1 Protein

Catalog No.: RP00202

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	6997	P13385

Tags

C-His

Synonyms

TDGF1;CR;CRGF;CRIPTO

Background

Cripto/TDGF1 is a member of the epidermal growth factor (EGF)- Cripto, Frl-1, and Cryptic (CFC) family. TDGF1 is an extracellular, membrane-bound signaling protein that plays an essential role in embryonic development and tumor growth. It is overexpressed in many types of cancers and acts as a growth factor for tumors. Cripto mutants display defects in mesoderm induction and heart morphogenesis, similar to phenotypes seen in Nodal mutants. Cripto can also activate mitogen-activated protein kinase (MAPK) and Akt pathways independently of Nodal by directly binding to a membrane-associated heparan sulfate proteoglycan, glypican-1.

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Basic Information

Description

Recombinant Human Cripto-1/TDGF1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Leu31-Thr172) of human Cripto/TDGF1 (Accession #NP_003203) fused with a 6×His tag at the C-terminus.

Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Recombinant Human Cripto at 2 μg/mL (100 μL/well) can bind Recombinant Human ALK4 with a linear range of 0.15-151 ng/mL.

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.


Centrifuge the vial before opening.
Reconstitute to a concentration of
0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting
the protein. For long term storage, it is
recommended to add a carrier protein
or stabilizer (e.g. 0.1% BSA, 5% HSA,
10% FBS or 5% Trehalose), and aliquot
the reconstituted protein solution to
minimize free-thaw cycles.

Contact




www.abclonal.com

Validation Data



Immobilized Recombinant Human Cripto at 2 µg/mL (100µL/well) can bind Recombinant Human ALK4 with a linear range of 0.15-151 ng/mL.



Recombinant Human Cripto-1/TDGF1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.