

RP00353

Leader in Biomolecular Solutions for Life Science



Recombinant Human Cyr61/CCN1 Protein

Catalog No.: RP00353 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	3491	O00622

Tags

C-hFc

Synonyms

Cellular communication network factor 1; Cysteine-rich angiogenic inducer 61; Insulin-like growth factor-binding protein 10; IBP-10; IGF-binding protein 10; IGFBP-10; Cyr61; CCN1; GIG1

Product Information

Source	Purification
HEK293 cells	≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

The secreted protein encoded by this gene is growth factor-inducible and promotes the adhesion of endothelial cells. The encoded protein interacts with several integrins and with heparan sulfate proteoglycan. This protein also plays a role in cell proliferation, differentiation, angiogenesis, apoptosis, and extracellular matrix formation.

Basic Information

Description

Recombinant Human Cyr61/CCN1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Thr25-Asp381) of human Cyr61/CCN1 (Accession #O00622) fused with a hFc tag at the C-terminus.

Bio-Activity

Measured by its ability to mediate Balb/3T3 mouse embryonic fibroblast cell adhesion. Recombinant Human Cyr61/CCN1 Fc Chimera immobilized at 1.5 μg/mL (100 μL/well), will mediate >25% BALB/3T3 cell adhesion (3 x 10⁴ cell/well added).

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human Cyr61/CCN1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.