

RP00364

Leader in Biomolecular Solutions for Life Science



Recombinant Human TNFRSF1B/TNF-R2/CD120b Protein

Catalog No.: RP00364 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	7133	P20333-1

Tags

C-His

Synonyms

CD120b; Etanercept; p75TBPII;
p75TNFR; TNF RII; TNF-R2; TNF-R75;
TNFR80; TNFRSF1B; TNFR2; TBPII;
TNFBR; TNFR1B

Product Information

Source	Purification
--------	--------------

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.	≥ 95 % as determined by SDS-PAGE; ≥ 95 % as determined by HPLC.
---	---

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Background

Tumor Necrosis Factor Receptor II (TNF RII), also known as TNFRSF1B, p75/p80, and CD120b, is a type I transmembrane protein that belongs to the TNF receptor superfamily. It has a molecular weight of approximately 75 kDa. Receptor with high affinity for TNFSF2/TNF-alpha and approximately 5-fold lower affinity for homotrimeric TNFSF1/lymphotoxin-alpha. The TRAF1/TRAF2 complex recruits the apoptotic suppressors BIRC2 and BIRC3 to TNFRSF1B/TNFR2. This receptor mediates most of the metabolic effects of TNF-alpha. Isoform 2 blocks TNF-alpha-induced apoptosis, which suggests that it regulates TNF-alpha function by antagonizing its biological activity.

Basic Information

Description

Recombinant Human TNFRSF1B/TNF-R2/CD120b Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Leu23-Asp257) of Human TNFR2/CD120b/TNFRSF1B (Accession #P20333-1) fused with a C-His tag at the C-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.

Reconstitution


Centrifuge the vial before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact




www.abclonal.com


Validation Data




Recombinant Human TNFRSF1B/TNF-R2/CD120b Protein was determined by Tris-Bis PAGE under reducing conditions.



Immobilized Human TNFR2, His Tag at 0.5 µg/mL (100 µL/Well). Dose response curve for Anti-TNFR2 Antibody, hFc Tag with the EC₅₀ of 6.4ng/mL determined by ELISA.



Immobilized Human TNFR2, His Tag at 5 µg/mL (100 µL/Well). Dose response curve for Biotinylated Human TNF alpha, His Tag with the EC₅₀ of 14.6ng/mL determined by ELISA (QC Test).



The purity of Human TNFR2 is greater than 95% as determined by SEC-HPLC.