

RP00436

Leader in Biomolecular Solutions for Life Science



Recombinant Human IL-15RA/CD215 Protein

Catalog No.: RP00436

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	3601	Q13261-1

Tags

C-His&Avi

Synonyms

IL-15 R alpha; CD215; IL-15RA;
MGC104179

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by SDS-PAGE; ≥ 95 % as determined by HPLC.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

Reconstitution

Background

Interleukin-15 receptor alpha (IL-15R alpha) is a high affinity IL-15 binding protein that is crucial for mediating IL-15 functions such as memory CD8 T cell proliferation and NK, NK/T cell, and intestinal intraepithelial lymphocyte development.

Basic Information

Description

Recombinant Human IL-15RA/CD215 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Ile31-Thr205) of Human IL-15RA/IL-15 R alpha/CD215 (Accession #Q13261-1) fused with a C-His&Avi tag at the C-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.


Centrifuge the vial before opening.
Reconstitute to a concentration of
0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting
the protein. For long term storage, it is
recommended to add a carrier protein
or stabilizer (e.g. 0.1% BSA, 5% HSA,
10% FBS or 5% Trehalose), and aliquot
the reconstituted protein solution to
minimize free-thaw cycles.

Contact




www.abclonal.com


Validation Data



The purity of Human IL-15RA is greater than 95% as determined by SEC-HPLC.



Recombinant Human IL-15RA/CD215 Protein was determined by Tris-Bis PAGE under reducing conditions.



Immobilized Human IL-15 at 1 µg/mL (100 µL/Well) on the plate. Dose response curve for Human IL-15RA, His Tag with the EC₅₀ of 22.1ng/mL determined by ELISA.