

RP00606

Leader in Biomolecular Solutions for Life Science



Recombinant Human Transcriptional Repressor Ctcf/CTCF Protein

Catalog No.: RP00606 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
<i>E. coli</i>	10664	P49711

Tags

No tag

Synonyms

CTCF; MRD21; CCCTC-binding factor;MRD21;CTCF

Product Information

Source

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Background

Transcriptional Repressor CTCF (CTCF) belongs to the CTCF Zinc-Finger Protein family. CTCF contains twelve C2H2-type zinc fingers and interacts with CHD8. CTCF is widely expressed in many tissues, and it is absent in primary spermatocytes. CTCF is involved in transcriptional regulation by binding to chromatin insulators and preventing interaction between promoter and nearby enhancers and silencers. CTCF plays an essential role in oocyte and preimplantation embryo development by activating or repressing transcription. In addition, CTCF is also indispensable in the epigenetic regulation and chromatin remodeling.

Basic Information

Description

Recombinant Human Transcriptional Repressor Ctcf/CTCF Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Met1-Ile154) of Human Transcriptional Repressor Ctcf/CTCF (Accession #P49711) fused with no tag.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.

Reconstitution


Centrifuge the tube before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human Transcriptional Repressor Ctcf/CTCF Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.