

RP01102

Leader in Biomolecular Solutions for Life Science



Recombinant Human IL-20RB Protein

Catalog No.: RP01102

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	53833	Q6UXL0

Tags

C-Fc

Synonyms

IL20RB;DIRS1;FNDC6;IL-20R2

Product Information

Source

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Supplied as a 0.22 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Background

IL20RB and IL20RA (MIM 605620) form a heterodimeric receptor for interleukin-20 (IL20; MIM 605619) (Blumberg et al., 2001 [PubMed 11163236]).

Basic Information

Description

Recombinant Human IL-20RB Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Asp30-Ala230) of human IL-20RB (Accession #Q6UXL0) fused with an Fc tag at the C-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.

Reconstitution


Centrifuge the tube before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human IL-20RB Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.