

RP01201

Leader in Biomolecular Solutions for Life Science



Recombinant Human ADAM8/CD156a Protein

Catalog No.: RP01201

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	101	P78325

Tags

C-His

Synonyms

ADAM8;CD156;CD156a;MS2

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of 50mM Tris, 10mM CaCl₂, 150mM NaCl, pH 7.5. Contact us for customized product form or formulation.

Background

Basic Information

Description

Recombinant Human ADAM8/CD156a Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Pro497) of human ADAM8/CD156a (Accession #NP_001100.3) fused with a 6×His tag at the C-terminus.

Bio-Activity

Measured by its ability to cleave a fluorogenic peptide substrate Mca-PLAQAV-Dpa-RSSSR-NH₂. The specific activity is >2 pmol/min/μg.

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.

Reconstitution


Centrifuge the vial before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human
ADAM8/CD156a Protein was
determined by SDS-PAGE under
reducing conditions with Coomassie
Blue.