

RP01335LQ

Leader in Biomolecular Solutions for Life Science



Recombinant Human TET2 Protein

Catalog No.: RP01335LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	54790	Q6N021-1

Tags

Please contact us for more information.

Synonyms

MDS; KIAA1546;TET2;KIAA1546

Product Information

Source	Purification
HEK293 cells	≥ 90 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Supplied as a 0.22 µm filtered solution of 50mM HEPES,300mM NaCl,1mM TCEP,10% Glycerol,pH 8.0.

Reconstitution

Background

Basic Information

Description

Recombinant Human TET2 Protein is produced by HEK293 cells expression system.

Bio-Activity

Measured by its binding ability in a functional ELISA. Immobilized Human TET2 Protein at 1 µg/mL (100 µL/well) can bind TET2 Rabbit pAb with a linear range of 0.976-232.5 ng/mL.

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -80°C
Avoid repeated freeze/thaw cycles.


Contact




www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human TET2 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.



Immobilized Recombinant Human TET2 Protein at 1μg/mL (100 μL/well) can bind TET2 Rabbit pAb with a linear range of 0.976-232.5ng/mL.