

RP01390

Leader in Biomolecular Solutions for Life Science



Recombinant Human Osteoactivin/GPNMB Protein

Catalog No.: RP01390

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	10457	Q14956-2

Tags

C-His

Synonyms

GPNMB;HGFIN;NMB

Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Endotoxin

<0.1EU/μg

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

Transmembrane glycoprotein NMB (GPNMB) is also known as Transmembrane glycoprotein HGFIN, DC-HIL and Osteoactivin (OA), which belongs to the PMEL/NMB family. GPNMB contains one PKD domain. GPNMB protein acts as a downstream mediator of BMP-2 effects on osteoblast differentiation and function. GPNMB participates in bone mineralization and functions as a negative regulator of inflammation in macrophages. Osteoactivin is expressed at high levels in normal and inflammatory liver macrophages suggesting a significant role in acute liver injury. The early-phase upregulation of Osteoactivin expression in the tubular epithelium in response to renal injury might play a role in triggering renal interstitial fibrosis via activation of matrix metalloproteinase expression and collagen remodeling in rats. Osteoactivin is a protein that is expressed in aggressive human breast cancers and is capable of promoting breast cancer metastasis to bone.

Basic Information

Description

Recombinant Human Osteoactivin/GPNMB Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Pro474) of human Osteoactivin/GPNMB (Accession #NP_002501.1) fused with a 6×His tag at the C-terminus.

Bio-Activity


Measured by its binding ability in a functional ELISA. Immobilized Human GPNMB Protein at 1 μg/mL (100 μL/well) can bind GPNMB Rabbit pAb with a linear range of 0.976-116.5 ng/mL.

Storage


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human Osteoactivin/GPNMB Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 70-100kDa.



Immobilized recombinant Human GPNMB Protein at 1µg/mL (100 µL/well) can bind GPNMB Rabbit pAb with a linear range of 0.976-116.5ng/mL.