

RP01456

Leader in Biomolecular Solutions for Life Science



Recombinant Human GSTA1 Protein

Catalog No.: RP01456 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	2938	P08263

Tags

C-His

Synonyms

GSTA1;GST-epsilon;GST2;GSTA1-1;GTH1

Product Information

Source	Purification
HEK293 cells	≥ 90 % as determined by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4, 5% trehalose, mannitol and 0.01% Tween80.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

GSTA1 (Glutathione S-Transferase Alpha 1) is a Protein Coding gene. This gene encodes a member of a family of enzymes that function to add glutathione to target electrophilic compounds. Glutathione S-transferases (GSTs) are involved in the detoxification of carcinogens and may be linked to carcinogenesis. As a vital component of GSTs, GSTA1 plays an important role in carcinogenesis. GSTA1 expression may be a target molecule in the early diagnosis and treatment of lung cancer. Human colonic adenocarcinoma (Caco-2) cells in culture undergo spontaneous differentiation into mature enterocytes in association with progressive increases in expression of glutathione S-transferase alpha-1 (GSTA1). GSTA1 levels may play a role in modulating enterocyte proliferation but do not influence differentiation or apoptosis. GSTA1 may play a key role during pregnancy.

Basic Information

Description

Recombinant Human GSTA1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Phe222) of human GSTA1 (Accession #NP_665683.1) fused with a 6×His tag at the C-terminus.

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human GSTA1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.