

RP01580

Leader in Biomolecular Solutions for Life Science



Recombinant Mouse SIRP-alpha/CD172a Protein

Catalog No.: RP01580

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	19261	P97797

Tags

C-His

Synonyms

CD172a; BIT; MFR; MYD1; MYD-1; P84; PTPNS1; SHP substrate 1; SHPS1; SHPS-1; SHPS1CD172A; SIRP alpha; SIRPA; Sirp-alpha-1; SIRPalph2; Sirp-alpha-2; Sirp-alpha-3; SIRPA

Product Information

Source	Purification
HEK293 cells	> 95% by SDS-PAGE.

Endotoxin

<0.1EU/μg

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Background

Signal regulatory protein α (SIRPα) is a regulatory membrane glycoprotein from SIRP family expressed mainly by myeloid cells and also by stem cells or neurons. SIRPα acts as inhibitory receptor and interacts with a broadly expressed transmembrane protein CD47 also called the "don't eat me" signal. Cancer cells highly expressed CD47 that activate SIRP α and inhibit macrophage-mediated destruction.

Basic Information

Description

Recombinant Mouse SIRP-alpha/CD172a Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Lys32-Asn373) of mouse SIRP-alpha/CD172a (Accession #NP_001277948.1) fused with a 6×His tag at the C-terminus.

Bio-Activity

Storage


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Mouse SIRP-
alpha/CD172a Protein was
determined by SDS-PAGE with
Coomassie Blue, showing a band at
60-70 kDa.