

RP01586LQ

Leader in Biomolecular Solutions for Life Science



# Recombinant Human Lysyl hydroxylase 2/PLOD2 Protein

Catalog No.: RP01586LQ **Recombinant**

## Sequence Information

| Species      | Gene ID | Swiss Prot |
|--------------|---------|------------|
| HEK293 cells | 5352    | O00469-1   |

### Tags

C-His

### Synonyms

PLOD2;BRKS2;LH2;TLH

## Product Information

| Source       | Purification                      |
|--------------|-----------------------------------|
| HEK293 cells | ≥ 95 % as determined by SDS-PAGE. |

### Endotoxin

< 0.1 EU/μg of the protein by LAL method.

### Formulation

Supplied as a 0.22 μm filtered solution in PBS, pH 7.4.

### Reconstitution

## Background

## Basic Information

### Description

Recombinant Human Lysyl hydroxylase 2/PLOD2 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (ALa26-Pro737) of human Lysyl hydroxylase 2/PLOD2 (Accession #NP\_000926.2) fused with a 6×His tag at the C-terminus.

### Bio-Activity

### Shipping

Lyophilized proteins are shipped at ambient temperature.  
Liquid proteins are shipped at frozen temperature with dry ice.  
Upon receipt, store it immediately at the temperature recommended below.

### Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

### Storage

Store at -80°C  
Avoid repeated freeze/thaw cycles.

## Contact




[www.abclonal.com](http://www.abclonal.com)

\* For your safety and health, please wear a lab coat and disposable gloves when handling.

## Validation Data

---



Recombinant Human Lysyl  
hydroxylase 2/PLOD2 Protein was  
determined by SDS-PAGE under  
reducing conditions with Coomassie  
Blue.