

RP01596S

Leader in Biomolecular Solutions for Life Science



Recombinant Mouse IL-13 Protein

Catalog No.: RP01596S

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	16163	P20109

Tags

NO-Tag

Synonyms

Interleukin-13; IL-13; T-cell activation protein P600 , IL13; IL-13

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 90% as determined by reducing SDS-PAGE.

Endotoxin

< 0.01 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of

Background

Mouse interleukin 13 (mIL-13) is a pleiotropic cytokine produced by activated Th2 cells. IL-13 induces B cell proliferation and immunoglobulin production. It contains a four helical bundle with two internal disulfide bonds. Mouse IL13 shares 58% sequence identity with human protein and exhibits cross-species activity. IL13 signals via receptor IL13R (type2, IL4R) and activates STAT-6. IL13 initially binds IL-13Rα1 with low affinity and triggers association of IL4R α , generating a high affinity heterodimeric receptor IL13R and eliciting downstream signals. IL13 also binds IL-13Rα2 with high affinity, which plays a role in a negative feedback system of IL13 signaling. IL13 is an important mediator of allergic inflammation and disease.

Basic Information

Description

Recombinant Mouse IL-13 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Gly21-Phe131) of Mouse IL-13 (Accession #P20109) fused with a His at the N-terminus.

Bio-Activity

Measured in a cell proliferation assay using TF-1 Human erythroleukemic cells. The ED₅₀ for this effect is 5.2-20.8 ng/mL, corresponding to a specific activity of 4.81×10⁴~1.44×10⁵ units/mg.

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.


0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact




www.abclonal.com

Validation Data



Measured in a cell proliferation assay using TF-1 Human erythroleukemic cells. The ED_{50} for this effect is 5.2-20.8 ng/mL, corresponding to a specific activity of $4.81 \times 10^4 \sim 1.44 \times 10^5$ units/mg.



Recombinant Mouse IL-13 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.