

RP01795

Leader in Biomolecular Solutions for Life Science



Recombinant Human S100-A7 Protein

Catalog No.: RP01795 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	6278	P31151

Tags

C-hFc

Synonyms

PSOR1; S100A7c; S100-A7; S100A7

Product Information

Source	Purification
HEK293 cells	> 97% by SDS-PAGE.

Endotoxin

< 0.01EU/μg

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

Protein S100-A7, also known as S100 calcium-binding protein A7, Psoriasin, S100A7, and PSOR1, is a secreted protein which belongs to the S-100 family. S100A7 was first isolated from skin involved by psoriasis, which can be induced in cultured squamous epithelial cells. S100A7 is expressed by both normal cultured and malignant keratinocytes and malignant breast epithelial cells within ductal carcinoma in situ, suggesting an association with abnormal pathways of differentiation. S100A7 plays a role in the pathogenesis of inflammatory skin disease, as a chemotactic factor for hematopoietic cells. It also plays a role in early stages of breast tumor progression in association with the development of the invasive phenotype.

Basic Information

Description


Recombinant Human S100-A7 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Ser2-Gln101) of human S100-A7 (Accession #NP_002954.2) fused with a hFc tag at the C-terminus.

Bio-Activity

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

Validation Data



Recombinant Human S100-A7 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 35-45 kDa.