

RP02006

Leader in Biomolecular Solutions for Life Science



Recombinant Human CD3 delta Protein

Catalog No.: RP02006

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	915	P04234

Tags

C-His

Synonyms

CD3D;CD3-DELTA;IMD19;T3D; IMD19; T3D

Product Information

Source

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by SDS-PAGE; ≥ 95 % as determined by HPLC.

Endotoxin

< 0.1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Background

T-cell surface glycoprotein CD3 delta chain, also known as CD3D, is a single-pass type I membrane protein. In immunology, the CD3 (cluster of differentiation 3) T cell co-receptor helps to activate both the cytotoxic T cell (CD8+ naive T cells) and also T helper cells (CD4+ naive T cells). It consists of a protein complex and is composed of four distinct chains. In mammals, the complex contains a CD3γ chain, a CD3δ chain, and two CD3ε chains.

Basic Information

Description

Recombinant Human CD3 delta Protein is produced by mammalian expression system. The target protein is expressed with sequence (Phe22-Ala105) of Human CD3 delta (Accession #P04234) fused with a 6xHis tag at the C-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C. Avoid repeated freeze/thaw cycles.

Reconstitution


Centrifuge the tube before opening.
Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact




www.abclonal.com

Validation Data



The purity of Human CD3 delta is greater than 95% as determined by SEC-HPLC.



Recombinant Human CD3 delta Protein was determined by Tris-Bis PAGE under reducing conditions.