

RP02082

Leader in Biomolecular Solutions for Life Science



Recombinant Human TNFRSF10A/DR4/TRAIL-R1/CD261 Protein

Catalog No.: RP02082

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	8797	O00220

Tags

C-His

Synonyms

TNFRSF10A;APO2;CD261;DR4;TRAILR-1;TRAILR1

Product Information

Source

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 90% as determined by SDS-PAGE.

Endotoxin

< 0.01 EU/μg of the protein by LAL method

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Background

Tumor necrosis factor receptor superfamily member 10A (TNFRSF10A) is also known as TNF-related apoptosis-inducing ligand receptor 1 (TRAIL-R1), Death receptor 4 (DR4), CD261 and APO2, which belongs to TNF superfamily. The expression of apoptosis-inducing TRAIL-R1 and TRAIL-R2 and of the decoy receptors TRAIL-R3 and TRAIL-R4 was systematically studied in all developmental stages of peripheral B cells isolated from the blood and secondary lymphoid organs. Expression of TRAIL-Rs is modulated along development, with highest levels observed in germinal center B cells.

Basic Information

Description

Recombinant Human TNFRSF10A/DR4/TRAIL-R1/CD261 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala109-Asn239) of Human TRAIL R1 / DR4 / TNFRSF10A fused with His tag at the C-terminal.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.


Centrifuge the tube before opening.
Reconstitute to a concentration of
0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting
the protein. For long term storage, it is
recommended to add a carrier protein
or stabilizer (e.g. 0.1% BSA, 5% HSA,
10% FBS or 5% Trehalose), and aliquot
the reconstituted protein solution to
minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human
TNFRSF10A/DR4/TRAIL-R1/CD261
Protein was determined by SDS-
PAGE under reducing conditions with
Coomassie Blue.