

RP02166

Leader in Biomolecular Solutions for Life Science



Recombinant Human Secretagogin/SCGN Protein

Catalog No.: RP02166 **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	10590	O76038

Tags

C-His

Synonyms

CALBL; DJ501N12.8; SECRET; SEGN; setagin; SCGN; DJ501N12.8; SECRET; SEG N; setagin

Product Information

Source	Purification
HEK293 cells	≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize freeze-thaw cycles.

Contact



www.abclonal.com

Background

Secretagogin (SCGN) is a secreted calcium-binding protein that is found in the cytoplasm; a small proportion is associated with secretory granules and membrane fractions. SCGN contains six EF-hand domains, related to calbindin D-28K and calretinin. SCGN is thought to be involved in KCL-stimulated calcium flux and cell proliferation. SCGN can be detected in human serum after ischemic neuronal damage. SCGN may function to negatively control growth and differentiation rates and, thus, indirectly inhibit cell replication.

Basic Information

Description

Recombinant Human Secretagogin/SCGN Protein is produced by Mammalian expression system. The target protein is expressed with sequence (Met1-Pro276) of human SCGN (Accession #O76038) fused with a 6xHis tag at the C- terminus.

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage


Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human Secretagogin/SCGN Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.