

RP02189

Leader in Biomolecular Solutions for Life Science



www.abclonal.com

Recombinant Human Galectin-7/LGALS7 Protein

Catalog No.: RP02189 **Recombinant**

Sequence Information

| Species | Gene ID | Swiss Prot |
|----------------|---------|------------|
| <I>E. coli</I> | 3963 | P47929 |

Tags

No tag

Synonyms

GAL7;LGALS7A;LGALS7

Product Information

| Source | Purification |
|----------------|-----------------------------------|
| <I>E. coli</I> | ≥ 95 % as determined by SDS-PAGE. |

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM EDTA, 5% Trehalose, pH8.0. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Background

The Galectin family of proteins, with specificity for N-acetylglucosamine containing glycoproteins, consists of beta-galactoside binding lectins containing homologous carbohydrate recognition domains (CRDs). They also possess hemagglutination activity, which is attributable to their bivalent carbohydrate binding properties. Galectins are active both intracellularly and extracellularly. Although they are localized primarily in the cytoplasm and lack a classical signal peptide; they can be secreted by one or more as yet unidentified non-classical secretory pathways. They have diverse effects on many cellular functions including adhesion, migration, polarity, chemotaxis, proliferation, apoptosis, and differentiation. Galectins may play a key role in many pathological states, including autoimmune diseases, allergic reactions, inflammation, tumor cell metastasis, atherosclerosis, and diabetic complications.

Basic Information

Description

Recombinant Human Galectin-7/LGALS7 Protein is produced by E.coli expression system. The target protein is expressed with sequence (Met1-Phe136) of human LGALS7 (Accession #P47929).

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human
Galectin-7/LGALS7 Protein was
determined by SDS-PAGE under
reducing conditions with Coomassie
Blue.