

RP02394B

Leader in Biomolecular Solutions for Life Science



Biotinylated Recombinant Human IL-13 Protein

Catalog No.: RP02394B

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	3596	AAK53823(P35225)

Tags

N-His&Avi

Synonyms

IL13;IL-13;P600

Product Information

Source

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by Tris-Bis PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the tube before opening.

Background

Interleukin-13 (IL-13) is a monomeric 17 kDa immunoregulatory cytokine that plays a key role in the pathogenesis of allergy, cancer, and tissue fibrosis. It is secreted by several helper T cell subsets, NK cells, mast cells, eosinophils, basophils, and visceral smooth muscle cells. Inhibits inflammatory cytokine production. Synergizes with IL2 in regulating interferon-gamma synthesis. May be critical in regulating inflammatory and immune responses.

Basic Information

Description

Biotinylated Recombinant Human IL-13 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gly21-Asn132) of Human IL-13 fused with a His tag and Avi at the N-terminal.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.



Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Immobilized Biotinylated Human IL-13, His Tag at 2µg/ml (100µl/well) on the streptavidin (5µg/ml) precoated plate. Dose response curve for Human IL-13Ra2, hFc Tag with the EC₅₀ of 1.4µg/ml determined by ELISA.

0