

RP02418

Leader in Biomolecular Solutions for Life Science



Recombinant Human IL-27RA/TCCR Protein

Catalog No.: RP02418

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	9466	Q6UWB1

Tags

C-hFc

Synonyms

IL27RA;CRL1;IL-27RA;IL27R;TCCR;WSX1;zcytor1

Product Information

Source

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 95 % as determined by Tris-Bis PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the tube before opening.

Background

Basic Information

Description

Recombinant Human IL-27 Ra / TCCR Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gly34-Lys516) of Human IL-27 Ra / TCCR fused with a hFc tag at the C-terminal.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.


Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human IL-27RA/TCCR
Protein was determined by Tris-Bis
PAGE under reducing conditions.