

RP02520

Leader in Biomolecular Solutions for Life Science



# Recombinant Mouse IL-17A/CTLA-8 Protein

Catalog No.: RP02520

Recombinant

1 Publications

## Sequence Information

Species	Gene ID	Swiss Prot
<i>E. coli</i>	16171	Q62386

### Tags

C-His

### Synonyms

IL17;Ctla8;IL-17;Ctla-8;IL-17A;IL17A;il17a

## Product Information

### Source

*E. coli*

### Purification

> 98% by SDS-PAGE.

### Endotoxin

<0.1EU/μg

### Formulation

Lyophilized from a 0.22 μm filtered solution of 20 mM sodium citrate, 0.2 M NaCl, pH 4.5.

### Reconstitution

Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water.

## Background

Interleukin-17A is a protein that in humans is encoded by the IL17A gene. The protein encoded by this gene is a proinflammatory cytokine produced by activated T cells. This cytokine regulates the activities of NF-kappaB and mitogen-activated protein kinases. This cytokine can stimulate the expression of IL6 and cyclooxygenase-2 (PTGS2/COX-2), as well as enhance the production of nitric oxide (NO).

## Basic Information

### Description

Recombinant Mouse IL-17A Protein is produced by Escherichia coli expression system. The target protein is expressed with sequence of Mouse IL-17A fused with polyhistidine tag at the C-terminus

### Bio-Activity

Measured by its ability to inhibit IL-6 secretion in 3T3 cells. The ED<sub>50</sub> for this effect is <1 ng/mL. The specific activity of recombinant mouse IL-17A is > 1 x 10<sup>6</sup> IU/mg.

### Storage

Store the lyophilized protein at -20°C to -80°C for long term. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week. Avoid repeated freeze/thaw cycles.

## Contact



[www.abclonal.com](http://www.abclonal.com)