

RP02762B

Leader in Biomolecular Solutions for Life Science



Biotinylated Recombinant Human LILRA4/ILT-7/CD85g Protein

Catalog No.: RP02762B

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	23547	P59901

Tags

C-His&Avi

Synonyms

CD85g; ILT7; ILT-7; ILT7MGC129598; LILRA4; MGC129597

Product Information

Source

Purification

≥ 95 % as determined by Tris-Bis PAGE; ≥ 95 % as determined by HPLC.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4. Normally 8% trehalose is added as protectant before lyophilization. Contact us for customized product form or formulation.

Reconstitution

Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Background

LILRA4, also known as ILT7 and CD85g, is an approximately 60-70 kDa variably glycosylated transmembrane protein that regulates immune cell activation. Mature human LILRA4 consists of a 423 amino acid (aa) extracellular domain (ECD) with four immunoglobulin-like domains, a 21 aa transmembrane segment, and a 32 aa cytoplasmic domain. LILRA4 function coreceptor to limit the innate immune responses to viral infections; signaling occurs via FCER1G (PubMed:16735691, PubMed:19564354). Down-regulates the production of IFNA1, IFNA2, IFNA4, IFNB1 and TNF by plasmacytoid dendritic cells that have been exposed to influenza virus or cytidine-phosphate-guanosine (CpG) dinucleotides, indicating it functions as negative regulator of TLR7 and TLR9 signaling cascades.

Basic Information

Description

Biotinylated Recombinant Human LILRA4/ILT-7/CD85g Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Glu24-Asn446) of Human LILRA4/ILT-7/CD85g (Accession #) fused with His tag and Avi tag at the C-terminal.

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt.

After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.


Avoid repeated freeze/thaw cycles.

Contact




www.abclonal.com

Validation Data



Biotinylated Recombinant Human LILRA4/ILT-7/CD85g Protein was determined by Tris-Bis PAGE under reducing conditions.



The purity of Biotinylated Human LILRA4 is greater than 95% as determined by SEC-HPLC.