

RP02791LQ

Leader in Biomolecular Solutions for Life Science



Recombinant Human Lactate Dehydrogenase A/LDHA Protein

Catalog No.: RP02791LQ

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
<i>E. coli</i>	3939	P00338

Tags

N-His

Synonyms

LDHM; GSD11; PIG19; HEL-S-133P;
LDHA

Product Information

Source	Purification
<i>E. coli</i>	≥ 95 % as determined by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Supplied as a 0.22 μm filtered solution in 20mM Tris, 250mM NaCl, 10% Glycerol, pH8.0.

Reconstitution

Background

L-Lactate Dehydrogenase A Chain (LDHA) is an enzyme that catalyzes the conversion of L-lactate and NAD⁺ to pyruvate and NADH in the final step of anaerobic glycolysis. LDHA contains an N-terminal coenzyme binding region, a central catalytic site, and at least nine utilized Lys acetylation and two Tyr phosphorylation sites. LDHA belongs to the lactate dehydrogenase family, expressed predominantly in muscle tissue. LDHA mutations have been linked to exertional myoglobinuria.

Basic Information

Description

Recombinant Human Lactate Dehydrogenase A/LDHA Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Met1-Phe332) of human Lactate Dehydrogenase A/LDHA (Accession #NP_005557.1) fused with a 6×His tag at the N-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -80°C
Avoid repeated freeze/thaw cycles.


Contact



www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human Lactate Dehydrogenase A/LDHA Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.