

RP02799

Leader in Biomolecular Solutions for Life Science



Recombinant Human CTHRC1 Protein

Catalog No.: RP02799

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
HEK293 cells	115908	Q96CG8

Tags

C-His

Synonyms

CTHRC1; UNQ762/PRO1550; CTHRC1

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 90% as determined by SDS-PAGE.

Endotoxin

< 1 EU/μg of the protein by LAL method.

Formulation

Lyophilized from a 0.22 μm filtered solution of PBS, pH 7.4.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of

Background

Collagen triple helix repeat-containing protein 1, also known as Protein NMTC1, and CTHRC1, is a secreted protein that is glycosylated and highly conserved from lower chordates to mammals. CTHRC1 expression was not detectable in normal arteries. However, it is transiently expressed in the arterial wall in response to injury where it may contribute to vascular remodeling by limiting collagen matrix deposition and promoting cell migration. A short collagen motif with 12 Gly-X-Y repeats appears to be responsible for trimerization of the CTHRC1 protein and this renders the molecule susceptible to cleavage by collagenase. CTHRC1 overexpression caused a dramatic reduction in collagen type I mRNA and protein levels. Currently available data indicate that Cthrc1 expression in vascular cells regulates transforming growth factor beta responsiveness, thereby impacting transforming growth factor beta target genes, including collagens. Additionally, CTHRC1 increases bone mass as a positive regulator of osteoblastic bone formation and offers an anabolic approach for the treatment of osteoporosis.

Basic Information

Description

Recombinant Human CTHRC1 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Ser31-Lys243) of human CTHRC1 (Accession #NP_612464.1) fused with a His tag at the C-terminus.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C. Avoid repeated freeze/thaw cycles.


0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human CTHRC1
Protein was determined by SDS-
PAGE under reducing conditions with
Coomassie Blue.