

RP02813

Leader in Biomolecular Solutions for Life Science



Recombinant Human VHL Protein

Catalog No.: RP02813

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
<i>E. coli</i>	7428	P40337

Tags

N-His&C-Avi

Synonyms

RCA1; VHL1; pVHL; HRCA1

Product Information

Source

Centrifuge the vial before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Purification

≥ 75 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Lyophilized from a 0.22 µm filtered solution of 50 mM Tris, pH 7.5.

Reconstitution

Centrifuge the vial before opening. Reconstitute to a concentration of

Background

Basic Information

Description

Recombinant Human VHL Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Met1-Asp213) of human VHL (Accession #NP_000542.1) fused with a 6×His tag at the N-terminus and an AVI tag at the C-terminus..

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature. Liquid proteins are shipped at frozen temperature with dry ice. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -20°C
Avoid repeated freeze/thaw cycles.


0.1-0.5 mg/mL in sterile distilled water.
Avoid vortex or vigorously pipetting the protein. For long term storage, it is recommended to add a carrier protein or stabilizer (e.g. 0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose), and aliquot the reconstituted protein solution to minimize free-thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant Human VHL Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.