

RP03165LQ

Leader in Biomolecular Solutions for Life Science



Recombinant Human PCBP1 Protein

Catalog No.: RP03165LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
<i>E. coli</i>	5093	Q15365

Tags

N-His

Synonyms

HNRPX; HNRPE1; hnRNP-X; HEL-S-85; hnRNP-E1;PCBP1

Product Information

Source

E. coli

Purification

≥ 85 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Supplied as sterile 50mM Tris, 500mM NaCl, 20% glycerol, pH 8.0

Reconstitution

Contact



www.abclonal.com

Background

Poly(rC)-binding protein 1, also known as Heterogeneous nuclear ribonucleoprotein E1, Alpha-CP1, Nucleic acid-binding protein SUB2.3 and PCBP1, is a family member of heterogeneous nuclear ribonucleoproteins (hnRNPs) that belong to RNA-binding proteins and bear three KH domains. PCBP1 is loosely bound in the nucleus. It may shuttle between the nucleus and the cytoplasm. It is abundantly expressed in skeletal muscle, thymus and peripheral blood leukocytes while a lower expression is observed in prostate, spleen, testis, ovary, small intestine, heart, liver, adrenal and thyroid glands. PCBP1 is widely expressed in many human tissues and involved in the regulation of transcription, transportation process, and function of RNA molecules. PCBP1 plays a pivotal role in post-transcriptional regulation for RNA metabolism and RNA function in gene expression. PCBP1 acts as a negative regulator of CD44 variants splicing in HepG2 cells, and loss of PCBP1 in human hepatic tumor contributes to the formation of a metastatic phenotype.

Basic Information

Description

Recombinant Human PCBP1 Protein is produced by *E. coli* expression system. The target protein is expressed with sequence (Asp 2-Cys 356) of human PCBP1 fused with a His tag at the N-terminus.

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes


For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human PCBP1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.