

RP03331LQ

Leader in Biomolecular Solutions for Life Science



Recombinant Human FGFR-2/KGFR/CD332 (V564F) Kinase

Catalog No.: RP03331LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
Baculovirus- Insect Cells	2263	P21802

Tags

N-His-StrepII

Synonyms

FGFR2; BEK; KGFR; KSAM; CD332;
Keratinocyte growth factor receptor;
Fibroblast growth factor receptor 2

Product Information

Source

Purification

≥ 90 % as
determined by
SDS-PAGE; ≥ 90 %
as determined by
HPLC.

Endotoxin

< 1 EU/μg of the protein by LAL
method.

Formulation

Supplied as a 0.22 μm filtered solution
in 50 mM Tris-HCl, 200 mM NaCl, 20%
glycerol, 1 mM DTT. (pH 7.5). Contact
us for customized product form or
formulation.

Reconstitution

Please use running water to thaw it
quickly.

Contact



www.abclonal.com

Background

FGFR2 also known as CD332 is a member of the fibroblast growth factor receptor (FGFR) family, which includes FGFR1-4, and FGFR1. FGFR1-4 are cell surface membrane receptors that possess tyrosine kinase activity. FGFR2 has two naturally occurring isoforms, FGFR2IIIb and FGFR2IIIc, created by splicing of the third immunoglobulin-like domain. And FGFR2 has important roles in embryonic development and tissue repair, especially bone and blood vessels. Like the other members of the fibroblast growth factor receptor family, these receptors signal by binding to their ligand and dimerisation (pairing of receptors), which causes the tyrosine kinase domains to initiate a cascade of intracellular signals.

Basic Information

Description

Recombinant Human FGFR-2/KGFR/CD332 (V564F) Protein is produced by Baculovirus-Insect Cells expression system. The target protein is expressed with sequence (Asp456-Glu768 (V564F)) of Human FGFR2 (Accession #P21802) fused with a N-His-StrepII tag.

Bio-Activity

The activity of FGFR2 is based on the MSA technology, and the content and ratio of the substrate and the product are directly separated and detected in real time and dynamically by the different migration rates of the substrate and the product after the enzymatic reaction.

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage


Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.

Aliquots below 10 μL are not advisable. Product must not be stored in diluted solutions. Avoid repeated freeze-thaw cycles.


Avoid repeated freeze/thaw cycles.

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human FGFR-2/KGFR/CD332 (V564F) Kinase was determined by SDS-PAGE under reducing (R) and non-reducing (NR) conditions.



The activity of FGFR2 is based on the MSA technology, and the content and ratio of the substrate and the product are directly separated and detected in real time and dynamically by the different migration rates of the substrate and the product after the enzymatic reaction.