

RP03348LQ

Leader in Biomolecular Solutions for Life Science



www.abclonal.com

# Recombinant Human PKAC gamma/PRKACG Kinase

Catalog No.: RP03348LQ **Recombinant**

## Sequence Information

Species	Gene ID	Swiss Prot
Baculovirus- Insect Cells	5568	P22612

### Tags

No tag

### Synonyms

PRKACG; PKA C-gamma; cAMP-dependent protein kinase catalytic subunit gamma

## Product Information

### Source

### Purification

≥ 90 % as determined by SDS-PAGE; ≥ 90 % as determined by HPLC.

### Endotoxin

< 1 EU/μg of the protein by LAL method.

### Formulation

Supplied as a 0.22 μm filtered solution in 50 mM Tris-HCl, 200 mM NaCl, 20% glycerol, 1 mM DTT. (pH 7.5). Contact us for customized product form or formulation.

### Reconstitution

Please use running water to thaw it quickly.

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Background

cAMP-dependent protein kinase catalytic subunit gamma is an enzyme that in humans is encoded by the PRKACG gene. Cyclic AMP-dependent protein kinase (PKA) consists of two catalytic subunits and a regulatory subunit dimer. This gene encodes the gamma form of its catalytic subunit. The gene is intronless and is thought to be a retrotransposon derived from the gene for the alpha form of the PKA catalytic subunit. And PRKACG has been shown to interact with Ryanodine receptor 2.

## Basic Information

### Description

Recombinant Human PKAC gamma/PRKACG Protein is produced by Baculovirus-Insect Cells expression system. The target protein is expressed with sequence (Gly2-Phe351) of Human PRKACG (Accession #P22612) fused with No tag.

### Bio-Activity

The activity of PRKACG is based on the MSA technology, and the content and ratio of the substrate and the product are directly separated and detected in real time and dynamically by the different migration rates of the substrate and the product after the enzymatic reaction.

### Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

### Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

### Storage


Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature.

Aliquots below 10 μL are not advisable. Product must not be stored in diluted solutions. Avoid repeated freeze-thaw cycles.


Avoid repeated freeze/thaw cycles.

\* For your safety and health, please wear a lab coat and disposable gloves when handling.

## Validation Data



Recombinant Human PKAC gamma/PRKACG Kinase was determined by SDS-PAGE under reducing (R) and non-reducing (NR) conditions.



The activity of PRKACG is based on the MSA technology, and the content and ratio of the substrate and the product are directly separated and detected in real time and dynamically by the different migration rates of the substrate and the product after the enzymatic reaction.