

RP10076LQ

Leader in Biomolecular Solutions for Life Science



Recombinant Human Ubiquitin (K48 only) protein

Catalog No.: RP10076LQ **Recombinant**

Sequence Information

Species **Gene ID** **Swiss Prot**

<I>E.
coli</I>

Tags

His Tag

Synonyms

Product Information

Source

<I>E. coli</I>

Purification

≥ 95 % as
determined by
SDS-PAGE.

Endotoxin

Please contact us for more
information.

Formulation

20 mM Tris, 150 mM NaCl, 2 mM βME,
10% Glycerol

Reconstitution

Background

Ubiquitin (Ub) is a 76 amino acid protein widely expressed in the cytoplasmic and nucleus of cells. Ub is posttranslationally conjugated to proteins by the E1, E2, E3 protein ubiquitination cascade. Ub can be conjugated on proteins as monoUb or polyUb chains. Protein ubiquitination plays both proteolytic and nonproteolytic functions. Usually, polyubiquitinated proteins are targeted to the 26S proteasome for proteolysis. Typical concentration to support in vitro ubiquitination is 50-100 μM. In Ub(K48 only), all lysine residues but lysine 48 were substituted with arginines, thereby it only supports the formation of lysine 48-linked polyubiquitin chains.

Basic Information

Description

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.


Contact



www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human Ubiquitin (K48 only) protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.