

RP10082LQ

Leader in Biomolecular Solutions for Life Science



Recombinant Human Ubiquitin (L73P) protein

Catalog No.: RP10082LQ **Recombinant**

Sequence Information

Species **Gene ID** **Swiss Prot**
<I>E. coli</I>

Tags
His Tag

Synonyms

Product Information

Source **Purification**
<I>E. coli</I> ~95% by SDS-PAGE

Endotoxin

Formulation
20 mM Tris, 150 mM NaCl, 2 mM βME,
10% Glycerol

Reconstitution

Background

Ubiquitin (Ub) is a 76 amino acid protein widely expressed in the cytoplasmic and nucleus of cells. Ub is posttranslationally conjugated to proteins by the E1, E2, E3 protein ubiquitination cascade. Ub can be conjugated on proteins as monoUb or polyUb chains. Protein ubiquitination plays both proteolytic and nonproteolytic functions. Usually, polyubiquitinated proteins are targeted to the 26S proteasome for proteolysis. Typical concentration to support in vitro ubiquitination is 50-100 μM. Studies have shown that polyUb chains formed by the Ub(L73P) mutant are resistant to deubiquitination catalyzed by many deubiquitination enzymes.

Basic Information

Description

Bio-Activity

Storage


Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant human Ubiquitin (L73P) protein was determined by SDS-PAGE with Coomassie Blue.