

RP10100

Leader in Biomolecular Solutions for Life Science



Recombinant human Deubiquitinase DESI2 protein

Catalog No.: RP10100

Recombinant

Sequence Information

| Species | Gene ID | Swiss Prot |
|---------|---------|------------|
| Human | | |

Tags
GST Tag

Synonyms

Product Information

| Source | Purification |
|----------------|--------------|
| <I>E. coli</I> | |

Endotoxin

Formulation

Supplied in 21 mM Tris, 150 mM NaCl, 2 mM BME, 10% Glycerol. Contact us for customized product form or formulation.

Reconstitution

Background

DESI-2 (desumoylating isopeptidase 2, also called Permuted Papain fold Peptidases of DsRNA viruses and Eukaryotes (PPPDE1)) and DESI-1 (also called PPPDE2) contain a Ub-binding domain and a cysteine protease domain. They catalyze both desumolyating and deubiquitinating activities.

Basic Information

Description

Bio-Activity

Storage


This product is stable at $\leq -70^{\circ}\text{C}$ for up to 6 months from the date of receipt.
For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant human DESI-2 Protein with GST Tag was determined by SDS-PAGE with Coomassie Blue.