

RP10184LQ

Leader in Biomolecular Solutions for Life Science



www.abclonal.com

Recombinant Human Ubiquitin(L73P) protein

Catalog No.: RP10184LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
<I>E. coli</I>		P0CG47

Tags

No Tag

Synonyms

Ubiquitin

Product Information

Source	Purification
<I>E. coli</I>	≥ 95 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

20 mM Tris, 150 mM NaCl, 2 mM βME, 10% Glycerol

Reconstitution

Background

Ubiquitin (Ub) is a 76 amino acid protein widely expressed in the cytoplasmic and nucleus of cells. Ub is posttranslationally conjugated to proteins by the E1, E2, E3 protein ubiquitination cascade. Ub can be conjugated on proteins as monoUb or polyUb chains. Protein ubiquitination plays both proteolytic and nonproteolytic functions. Usually, polyubiquitinated proteins are targeted to the 26S proteasome for proteolysis. Typical concentration to support in vitro ubiquitination is 50-100 μM. Studies have shown that polyUb chains formed by the Ub(L73P) mutant are resistant to deubiquitination catalyzed by many deubiquitination enzymes.

Basic Information

Description

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.


Contact



www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant Human Ubiquitin(L73P) protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.