

RP10219K

Leader in Biomolecular Solutions for Life Science



Ubiquitin Linkage Screening Kit III

Catalog No.: RP10219K

Sequence Information

| Species | Gene ID | Swiss Prot |
|-------------|---------|------------|
| E. coli | | |
| recombinant | | |

Tags

Synonyms

Product Information

| Source | Purification |
|---------------------|-----------------------|
| E. coli recombinant | 95% based on SDS-PAGE |

Endotoxin

Formulation

20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol

Reconstitution

Background

Ubiquitin Linkage Screening Kit III is designed for two purposes: 1) to identify the ubiquitin lysine residue that is used to link the polyubiquitin chain conjugated on a protein substrate; 2) to determine whether a protein substrate is polyubiquitinated or conjugated with multiple monoubiquitin.

Basic Information

Description

Bio-Activity

Storage

Store at -70°C . This product is stable at $\leq -70^{\circ}\text{C}$ for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.