

RP10237K

Leader in Biomolecular Solutions for Life Science



www.abclonal.com

Deubiquitinating Enzyme (DUB) Identification Kit

Catalog No.: RP10237K

Sequence Information

Species	Gene ID	Swiss Prot
E. coli recombinant		

Tags

Synonyms

Product Information

Source	Purification
E. coli recombinant	≥ 95 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

Reconstitution

Background

This kit is designed to enrich deubiquitinating enzymes in crude or partially purified cell/tissue extracts, followed by identification using immunoblotting or mass spectrometry.

Basic Information

Description

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -70°C. This product is stable at $\leq -70^{\circ}\text{C}$ for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.