

RP10240K

Leader in Biomolecular Solutions for Life Science



www.abclonal.com

Caspase-3, 7 Activity Fluorometric Assay Kit

Catalog No.: RP10240K

Sequence Information

Species **Gene ID** **Swiss Prot**

Cell Lysates

Tags

Synonyms

Product Information

Source

Cell Lysates

Purification

Endotoxin

Please contact us for more information.

Formulation

Reconstitution

Background

Caspase-3, 7 Activity Fluorometric Assay Kit can be used for assaying caspase-3, 7 activities in cell/tissue extracts in a 96-well plate format. Each of the supplied fluorogenic substrate and the pan caspase inhibitor Z-VAD-FMK is sufficient for 250 x 100 µl reactions. The supplied TRAIL-treated cell lysate is used as a positive control, sufficient for 5 X 100 µl reactions.

Basic Information

Description

Bio-Activity

Shipping

The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -70°C. This product is stable at $\leq -70^{\circ}\text{C}$ for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.